

Programat P95


VDE

Operating Instructions

Schulungsunterlage
Course material

CE

IVOCLAR



Contents

	Page		Page
Supply packs	3	Removal of heating muffle	27
Important points to be noted	4	Assembly of heating muffle	29
Parts list	4	Furnace servicing	30
Pre-operating procedure	6	Vacuum pump servicing	30
Connections	7	Technical data	31
Switches, keypad and indicators	8	Summary of programs for operation in Celsius mode	32
Explanation of terms	10	After-sales service companies for assistance and repairs	33
Trial run	11		
Programs P01-P91 and P99	12		
Operation in either Celsius or Fahrenheit mode	14		
Controlling the CTE value	14		
Programming, change of program	15		
Important practical information	16		
Checking the furnace temperature (silver test)	17		
Recalibration	18		
Errors not indicated on the display and their correction	18		
Indicated errors and their correction	19		
Removal of control unit	23		
	24		
	25		
Change of heating muffle	26		

The Programat P95 is a fully automatic ceramic vacuum furnace with sophisticated microprocessor control. It was developed from the tried-and-tested Programat P80.

The temperature sensor has been greatly improved. The platinum tube now used provides optimum conduction and transfers the values to the microprocessor at a significantly higher speed.

In the area of software, practical experience has led to improvements which give the dental technician more individual options and more effective control over the ceramic material to be fired.

The Programat P95 is supplied in four module kits, each packed in polystyrene:

Supply packs

Pack 1:
Furnace head (K), cooling plate (A)

Pack 2:
Control unit (S)

Pack 3:
Furnace base (U), vacuum hose, spare fuses

Pack 4:
Accessories:
Set comprising tongs, firing trays K and G, temperature control set and program cards

In addition to the standard cream colour, Programat P95 is available in the following colours:


Colours	IVOCULAR colour name:	RAL number and name:
Standard colour	cream standard	RAL 1013 pearl white
Special colours:	salmon pink	RAL 5014 old rose
	aquamarine	RAL 5014 pigeon blue
	turquoise	RAL 6027 light green
	gray	RAL 7035 light gray
	white	RAL 9016 traffic-sign white

Important points to be noted

- Avoid positioning the furnace and pump in the immediate vicinity of radiators or other sources of heat.
- Ensure that the air vents in the rear panel of the furnace are unrestricted at all times.
- Install the vacuum pump in a well-ventilated place. Ensure that the apertures (68a, b, c) in the frame plate (68) are free and that no foreign matter can fall into the furnace base.
- Avoid resting any objects on the frame plate (68); rest objects only on the cooling plate (A).
- Ensure that the sealing ring (35) in the furnace head and the sealing rim (69) of the furnace base are kept clean and undamaged.
- Be careful not to touch the hot parts of the furnace during operation.
- Do not damage the blade contacts (33).
- Clean only with a dry or slightly moist cloth. Do not use solvents.
- For forwarding, use original packaging.


Parts list

- S-Control unit
- 1 0/1 switch with pilot lamp
 - 2 Keypad
 - 3 Display window
 - 4 Front panel
 - 5 Cover
 - 6 Tab
 - 7 Screws
 - 8 Tab
 - 9 Tab
 - 10 Battery
 - 11 Parallel-circuit board
 - 12 A-C cooling plate




K-Furnace head

- 30 Dome
- 30a Inside wall of hood
- 31 Mounting lug
- 31a Holes
- 32 Switch pin
- 33 Blade contacts
- 34 Protective caps
- 35 Sealing ring
- 36 Spring clip
- 37 Ring
- 38 Insulation cord
- 38a End
- 39 Stone lining segments
- 39a Joints
- 40 Heating muffle
- 40a,b End of wire
- 41a,b Terminal
- 42a,b Terminal screws



U-Furnace base

- 60 Casing
- 61 Front panel
- 61a Guides
- 61b Floor of base
- 61c Locking bolt
- 61d Pad
- 62 Vacuum gauge
- 63 Air vents
- 64 Mains lead
- 65 Pump power socket
- 66a Pump fuse
- 66b Mains fuse
- 66c Heating circuit fuse
- 67 Vacuum hose connection
- 68 Frame plate
- 68a,b,c Apertures in frame plate
- 69 Sealing rim
- 70 Stone lining insert
- 71 Sheathed thermocouple
- 72 Firing mount
- 73 Hinge pin
- 89 Pump plug
- 99 Vacuum hose



Pre-operating procedure

1. Unpacking

- Grip the individual parts through the recesses in the packaging and carefully remove from the packaging.
- Check for any damage in transit and clean off residues of polystyrene (keep the packaging).
- Check that the voltages indicated on the plates on the back of the furnace head, on the furnace base and on the control unit comply with the local mains voltage.

Important:


The sheathed thermocouple must be set perpendicular and must not be either damaged or bent.

2. Assembling the furnace base and its sections


- Position the furnace base (1)
- Remove the firing mount (72) from the packaging material and introduce it into the stone lining insert (73).
- Clean the sealing rim (69)
- Mount the cooling plate (A) in the three indentations in the frame plate (68).

3. Mounting the furnace head

- Blow out the muffle (40) and the surfaces of the stone lining (39) at moderately low pressure or clean carefully with a soft brush. Do not touch the heating element!
- Clean the sealing ring (35) of the furnace head, the heating muffle and the surfaces of the stone lining segments (39). Do not touch the heating element of the heating muffle!



- Set the firing pins (73) perpendicular to the furnace base.



- Hold the furnace head and fit through the apertures (31a) in the mounting lug (31) onto the firing pins (73).
- Keeping the furnace head level, push down in a parallel direction until the sealing ring (35) rests evenly on the sealing rim (69) of the furnace base.

Connections

1. Vacuum pump

- Programat vacuum pump:
Follow separate pump operating instructions.
 - Pumps from other manufacturers:
(for permitted specifications see page 32).
- a) Have the Programat pump plug connected to the vacuum pump mains lead by a qualified electrician.
 - b) Position the pump and insert the pump plug into the socket (65).
 - c) Push the vacuum hose into the hose connection (67) of the furnace and the vacuum pump connection.

Important:

Before operating the furnace, connect the vacuum pump to the base of the furnace!

2. Furnace

Connect the mains lead (64) to the mains supply.

Switches, keypad and indicators



1. 0/1 switch with green pilot lamp

When the 0/1 switch has been pressed, the green pilot lamp is lit, the buzzer sounds and various data appear briefly on the displays. After about six seconds, the effective values are displayed. The closed furnace heats to stand-by temperature B.

2. Selection keys

Key	To select:
P	Program
B	Stand-by temperature
T	Temperature increase
F	Firing temperature
S	Closing time
H	Holding time
L	Long-term cooling
V ₁	Vacuum on
V ₂	Vacuum off

When a selection key is pressed, a dot is lit at the end of the relevant display. This indicates data input or that the value can be changed.

3. Stand-by temperature

The stand-by temperature can be freely programmed from 100-700°C (212-1292 °F).

4. Holding time

Holding time is entered in tenths of a minute from 0.1 to 9.9 minutes.

Holding times of 10-99 minutes are entered in steps of one minute.

Example: holding time 1.3 - enter H1 H3 → display 1.3

Input values in seconds:

0.1 - 6 seconds

0.2 - 12 seconds

0.3 - 18 seconds

0.4 - 24 seconds

0.5 - 30 seconds etc.

Second holding time

Some ceramic compounds have to be processed with a second holding time during heating. Programs P70-P75 with an automatically set second holding time (see pages 12 and 31) have been provided for this purpose. **Holding time 50% with vacuum and 50% without vacuum.**

In program P76, the holding time is automatically divided:

50% = first half with vacuum, 50% = 2nd half without vacuum.

5. Data input keys

0 to 9: for the input of program number and desired values. For possible desired values, see pages 12 and 31.

6. Control keys

↑ : to open furnace

↓ : to close furnace

The control keys ↑ (open furnace) and ↓ (close furnace) cannot be operated while a program is running. If this is desired, however, first press the STOP key and then the control keys.

STOP: Press once to:

- stop a program which is running
- stop furnace head movement
- stop heating
- stop the buzzer

Press twice to:

- stop the vacuum

START: Press the green key to start a program.

7. Calibration keys

— and — calibration of furnace temperature (for instructions see page 18).

8. Indicators in the long displays

The displays show:

- Program number as selected with key R
- Desired value as selected with the relevant key (no data displayed for L, V₁, V₂ if data input = 0; 0.0 is displayed at H)
- Duration of program sequence is shown on the clock display once a program has been started. (Approximate time in minutes shown after the flashing dot. If L is pre-set below 500 °C [932 °F], only the dot flashes).

Important:

Flashing of a data indicator signals an error in input (for explanation see page 19).

9. Indicators in the small display

- Furnace temperature (effective T value)
- ▼ = program running
- Flashing = error indication (for explanation, see pages 19 to 22).

10. Vacuum indicator

- Indicator at end of scale = no vacuum
- Indicator in green sector = operating vacuum

11. Buzzer

The buzzer sounds to signal

- close the open furnace
- position or remove object on firing mount
- end of program

Time sequence: 10 seconds buzzing - approx. 5 minutes pause - approx. 5 minutes buzzing.

The acoustic signal can be changed using program P91. For description see page 13.

Explanation of terms

P = Program

B = Stand-by temperature

This is the temperature to which the furnace heats after it has been closed and switched on, but no program has been started.

IT = Temperature increase per minute

Increase in furnace temperature after the furnace has been closed and a program started.

T = Firing temperature

The furnace temperature held constant once the desired temperature increase I has taken place.

N = Duration of program in minutes with countdown.

S = Closing time (minutes)

Furnace closing time once the program has been started.

H = Holding time (minutes)

The period during which the furnace temperature is kept constant.

L = Long-term cooling

The furnace only opens when the temperature has fallen to the preset L value.

V₁ = Vacuum on

Development of vacuum (pump starts).

V₂ = Vacuum off

Release of vacuum (pump stops).



Trial run

1. Switch on

Press 0/1 switch and wait until the furnace reaches the stand-by temperature **B**.

2. Check operating mode

Set for operation in Celsius or Fahrenheit mode (for procedure, see page 14).

3. Enter trial program

First press program selection key, e.g. P09.

Program	°C mode	°F mode
B	400 (°C)	762 (°F)
φ	140 (°C/min.)	252 (°F/min.)
T	730 (°C)	1346 (°F)
S	04 (min.)	04 (min.)
H	09 (min.)	09 (min.)
L	650 (°C)	1234 (°F)
V ₁	510 (°C)	950 (°F)
V ₂	700 (°C)	1292 (°F)

4. Open furnace

Press key 1 (furnace open).

5. Run program

Press **START** after the buzzer sounds. The clock indicator shows the duration of the program sequence at any given time.


The program runs automatically (see diagram on the right).

6. Close furnace

The buzzer indicates the end of the program. Close the furnace with the key.

Note:

The outside of the furnace head becomes hot when the furnace is open!



Programs P01-P91 and P99

Programs P01 to P08

Standard programs with some data automatically set. The following desired values can be freely programmed within the limits of the T, S, H, L (in P07 and P08 also V₁ and V₂) listed on page 31 (Celsius mode).

Program P01

Heat treatment (oxidation) under vacuum.

Automatically set:

B = 403 °C (757 °F)

T = 140 °C / min. (252 °F / min.)

V₁ = As soon as the furnace is closed

V₂ = end of holding time

Program P02

As Program P01, but without vacuum.

Program P03

Vacuum firings (dentin, incisal)

Automatically set:

B = 403 °C (757 °F)

T = 60 °C / min. (108 °F / min.)

V = 580 °C (1076 °F)

V₂ = start of holding time

Program P04

Glaze firings, as P03, but without vacuum.

Program P05

Opaque firings (dentin, incisal)

Automatically set:

B = 403 °C (757 °F)

T = 80 °C / min. (144 °F / min.)

V₁ = 550 °C (1022 °F)

V₂ = start of holding time

Program P06

Glaze firings, but without vacuum

Automatically set:

B = 403 °C (757 °F)

T = 80 °C / min. (144 °F / min.)

Program P07

Dentin firings "V", with vacuum

Automatically set:

B = 403 °C (757 °F)

T = 60 °C / min. (108 °F / min.)

(V₁ and V₂ freely programmable)

Program P08

Opaque firings, with vacuum

Automatically set:

B = 403 °C (757 °F)

T = 80 °C / min. (144 °F / min.)

(V₁ and V₂ freely programmable)

Program P09 to P69

Freely programmable programs with normal opening of the furnace (1 minute).

For possible desired values, see page 31 (Celsius mode).

Programs P70 to P75

Special programs, each with a second, automatically set holding time and temperature.

For other possible desired values, see page 31 (Celsius mode).

Program	Automatically set holding temperature	Automatically set holding time
P70	575 °C (1067 °F)	2 minutes
P71	575 °C (1067 °F)	3 "
P72	600 °C (1112 °F)	2 "
P73	600 °C (1112 °F)	3 "
P74	625 °C (1157 °F)	2 "
P75	625 °C (1157 °F)	3 "

Program P76
Special program in which the first half of the holding time, which has to be entered specially, is carried out with vacuum and the second half without vacuum.

All values are freely programmable.
Example: holding time H3 carried out as 1.5 min. with vacuum, 1.5 minutes without vacuum.

Programs P77 to P87
Industrially programmable special programs with rapid opening of the furnace (20 seconds). For possible desired values, see page 31 (Celsius mode).

Program P88 to P90
Freely programmable 'night-time' programs with normal opening of the furnace (1 minute). After the 'night-time' program has been completed, the heating switches off without the buzzer sounding and the furnace closes automatically after reaching a temperature of approx. 80 °C (176 °F) and cools to room temperature.

In the event of a power failure during the night, the furnace does not continue heating, but remains at room temperature.

Program P91
With this program the acoustic buzzer signal can be individually set or completely switched off. Different buzzer signals can thus be generated if several devices are being operated simultaneously.

Set on-time with V₁ key (editing dot is lit) and V₄ or V₂ key from 0 to 9. The display in the V₁ sector corresponds to the on-time (20 msec.). Display 0 = buzzer completely switched off. Set interval with V₂ key (editing dot is lit) and V₁ or V₂ key from 0 to 9. The display in the V₂ sector corresponds to the off-time (200 msec.). Display 0 = buzzer sounds continuously. Standard setting: V₁ = 5 (100 msec. On) and V₂ = 1 (100 msec. Off).

Note:
Programming is only possible when the editing dot is lit.

Program P98
This program can be used to change the Program P95 from Celsius mode (°C) to Fahrenheit mode (°F) and vice-versa.

Program P99
Program for checking the furnace temperature with the 'silver test' and for 're-calibration'. For detailed description, see page 17.

Summary of programs for operation in Celsius mode

Possible desired values for operation in Celsius mode
(For list of programs see page 31.)

Operation in Celsius ($^{\circ}\text{C}$) or Fahrenheit ($^{\circ}\text{F}$) mode

Note:

Program P98 is not a working program. Error indication E100 and flashing of the program number are therefore of no significance.

1. Mode check

Enter P98:

- If a 'C' is lit above the V₁ key, the furnace is in Celsius mode.
- If an 'F' is lit above the V₁ key, the furnace is in Fahrenheit mode.

2. Changeover in Fahrenheit mode

Enter P98:

- C is lit above the V₁ key
- Switch off 0/1.
- Press V₁ and keep pressing while switching on 0/1.
- Release V₁ key approx. 3 seconds after switching on. As soon as the F lights up above the V₁ key, the furnace is in Fahrenheit mode.

3. Changeover from Fahrenheit to Celsius mode

Follow the procedure described under 2.

As soon the C lights up above the V₁ key, the furnace is in Celsius mode.

door is open (minus)	programs P9C to P97 (without long-term cooling L)	-0.5			remove object as soon as furnace is completely open or no later than when the first buzzer signal sounds
standard (0)	programs P01 to P99 (without long-term cooling L)	+0			remove object when the buzzer signals sound
upstarts (plus)	programs P01 to P99 (with long-term cooling L)	+0.5 L = 800°C (1472°F)	-1.0 L = 600°C (1122°F)	+1.5 L = 400°C (752°F)	remove object when the buzzer signals sound

Controlling the CTE value

The CTE (Coefficient of Thermal Expansion) value of the ceramic material can be controlled as follows:

1. Immediate removal of the object from the furnace after firing results in a decrease in CTE (minus).

2. Slow cooling of the object in the furnace after firing (long-term cooling) results in an increase in CTE (plus).

Programming, change of program

- Program cards are available for noting the program data.
- Program numbers should be given in two digits, i.e. P01 to P09, P10, P11. (Not P1, P2 etc.).
- As long as no program is running, data may be entered or modified by the following procedure:
- Press selection key (A dot is lit at the end of the relevant display).
- Enter the value.
- Important for the input of V_2 (vacuum off):
- If vacuum-off firing takes place during holding time H, the following input should be used with V_2 :
Celsius mode: $V_2 = T - 1^\circ\text{C}$
(e.g. $T = 1050^\circ\text{C}$, $V_2 = 1049^\circ\text{C}$)
Fahrenheit mode: $V_2 = T - 2^\circ\text{F}$
(e.g. $T = 1922^\circ\text{F}$, $V_2 = 1920^\circ\text{F}$)
(Vacuum is switched off at the start of holding time H).


- If vacuum-on firing takes place during the holding time, the following input should be used:
 $V_2 = T$ (e.g. $T = 1050^\circ\text{C}$, $V_2 = 1050^\circ\text{C}$
or $T = 1922^\circ\text{F}$, $V_2 = 1922^\circ\text{F}$)
(Vacuum is not switched off until the end of the holding time H).
 - If control is lost because of incorrect programming, the situation can be rectified by the following procedure:
 - Switch off D/I.
 - Press the STOP key and keep pressing while switching on again.
(The displays then show the desired values originally set by the manufacturer.)
 - Once the program has been completed, it is automatically stored.
- Program changeover is possible at any time while the program is running:
 - 1. To change from P01 to P02, P03 to P04, P05 to P06, P07 to P08 and vice-versa: press P and enter new program number (the program sequence is not interrupted).
 - 2. To change from one program to another:
 - Press STOP.
 - Press P and enter new program number.
 - Press START.
 - It is possible to change preselected data while a program is running only if the effective temperature has not yet reached the preselected I value:
 1. Desired values for S, H and L can be changed without interrupting the program sequence: Press the selection keys and enter the new value.
 2. To change preselected data for α , β , V_1 and V_2 : Press STOP.

- Press the appropriate selection key and enter the new value.
- Press START.
- To interrupt a vacuum program which is running and release the vacuum: Press STOP twice.

Important practical information

- Always keep the furnace closed between firings.
- Optimum results can be obtained with Ivoclar silicon nitride firing mounts.
- Objects which have to be pre-dried should be placed on the firing mount only after the buzzer has sounded.
- A power failure will interrupt a program which is running and cause Er17 to appear. Press STOP and re-start the program for the program to continue its sequence. (Any adverse effect on the object depends on how long the power failure lasts.)
- Check the furnace temperature by carrying out the 'silver test' (see page 17).
- After a control unit or furnace base has been changed, it is advisable to check the furnace temperature.
- Do not open the furnace head manually when the furnace is switched on. Error Er28 will otherwise be indicated. (For correction of error see page 22.)

Checking the furnace temperature (silver test)



The sheathed thermocouple may be subject to changes which affect the furnace temperature, depending on the mode and period of operation.


Material required:
(in the temperature checking set)

- Firing tray marked IVOCLAR-E (17 mm dia., white)
- Strip of silver, purity 99.99% (size 0.5 x 2 x 20 mm)

Procedure:

- a) The furnace must be at operating temperature (switched on for at least 30 minutes) and have a stand-by temperature of 400 °C (e.g. in P99).
- b) Insert silver strip into the IVOCLAR-E firing tray.
- c) Select P99 (silver test program).
- d) Press F key and place firing tray with silver strip in the centre of the firing mount (72).
- e) Press START. (If error indication Tr14 appears, the furnace temperature is still too high for the 'silver test'. The furnace closes automatically at the correct temperature and the program starts).


If the silver strip is melted at the end of the program, the furnace temperature is correctly calibrated. If not, recalibration is necessary.



Recalibration

A change in temperature of + 50 °C (= +90 °F) is possible in the Programat P95. Select program P99 to activate the calibration keys + and -.

- If the silver strip has not started to melt after the silver strip test, recalibrate using the + key.
- If the silver strip has melted down into a well after the silver test, recalibrate using the - key.



Every time a calibration key is pressed, the set temperature changes by 1 °C (1.8 °F). Experience has shown that a recalibration of 5 °C (9 °F) is appropriate, which means pressing the relevant calibration key five times. While the calibration keys are in use, the calibration value is shown in the small display in °C. Repeat the silver test until the silver strip has started to melt.

Errors not indicated on the display and their correction

- Green pilot lamp fails to light although the 0/1 switch is on:
 1. The indicators are not lit: check fuses (66a) and (66b); replace if defective.
 2. The indicators are lit:
Replace the defective pilot lamp. The furnace is nevertheless capable of operating.
- Incomplete or illogical values are displayed:
 - Switch off 0/1
 - Press STOP and keep pressing while switching 0/1 on again.If incomplete or illogical values continue to appear, the control unit is defective. Replace as described in the instructions (page 7 and 25)

- No vacuum, pump is running:
Check sealing rim (69) and sealing ring (35) and clean if necessary. Check that the vacuum hose (59) has no leaks and is firmly attached at both ends. If vacuum is still not observed, consult after-sales service.

- No vacuum, pump is not running:
Check fuses (66a and 66b). Replace if defective. If the pump still fails to operate, have it checked by an electrician. If the pump appears to be in order but still fails to operate, consult after-sales service.

Indicated errors Er00 to Er11 and their correction

Er00 to Er11 indicate incorrect input. Er. and the input to be corrected flash in the effective temperature indicator. (Program is interrupted).

Error	Cause of error	Correction
Er00	Incorrect P number 00 1 to 9 91 to 98	Press selection key P and enter correct program number*) Error indication in P98 indicates that this program is not a working program*)
Er01	T input below actual temperature (furnace opens)	Press selection key above which preselected value flashes and then enter correct value*.)
Er02	T input below R or above 1200°C (2192°F)	
Er03	S input incorrect 0/0, 0/0, 1/0, 2 minutes	
Er04	H input incorrect over 30 minutes	
Er05	r input: 0 to 29°C/min, or over 140°C/min, 0 to 53°F/min or above 252°F/min.]	

*) Take correct program number and correct desired value from the summary table on page 31 (Celsius mode).

Er06	B input 0 to 99 °C or above 700 °C (0 to 210 °F or above 1293 °F)	Press selection key above which preselected value flashes and then enter correct value.*
Er07	B input above desired T value	
Er08	L input above desired L value	
Er09	V ₁ input above desired V ₁ value	
Er10	V ₂ input above desired T value	
Er11	V ₁ input or V ₂ input absent	

Indicated errors Er12 to Er99 and their correction

Error indication Er_n appears in the effective temperature indicator.

Error	Cause of error	Correction
Er12	Battery in control unit is low and has to be replaced quickly with a new one.	Press STOP (Er12 goes over). Obtain new battery from the after-sales service unit and install it in accordance with the instructions on page 24. Until battery is changed, operate furnace as usual.
Er13	Overheating	1. Error in programming: most preselected T value lower than effective temperature of furnace.
Er14	Furnace chamber temperature too low, for Silver level	Furnace closes automatically at correct temperature and program runs.
Er17	Means failure when program has started	Press STOP, Er17 goes out. It is then possible to continue working.
Er20	Detect in furnace base or in control unit (furnace opens).	2
Er21	Detect in control unit (furnace opens)	1
Er22	Detect in furnace base or in control unit (furnace opens)?	2

1 → Switch off 0/1 and switch on again after a few seconds.
If Er_n re-appears, press STOP and repeat the procedure.
If Er_n continues to flash, replace control unit according to the instructions on pages 7 and 23. (State Er number when requesting repairs).


2 → Proceed as described under 1, but without replacing the control unit, and consult after-sales service.

Er23	Heating muffle defective (furnace does not heat)	Press STOP (Er 23 goes out). Work can proceed as normal. Obtain new heating muffle from after-sales service unit and install in accordance with instructions on pages 29 to 30.
Er24	Heating muffle defective (furnace does not heat)	Switch off. Obtain new heating muffle and install in accordance with instructions on pages 29 to 30.
Er25	Internal chamber temperature above 55°C (131°F)	1. Chamber temperature too high: cool furnace base with fan 2. Initialization, see page 25.
Er26	Loss of stove status	1. Initialization, see page 25.
Er27	Defect in control unit (furnace opens)	1
Er28	Heating element interrupted; furnace was opened manually or is defective	2. If furnace was not opened manually, check heating circuit fuse (Er5).
Er29	Defect in control unit	1. Initialization, see page 25.
Er30	Defect in control unit	1
Er32	Battery low	Replace battery and initialize see page 25.
Er60 to Er99	Defect in control unit	1

1 → Switch off 0/1 and switch on again after a few seconds.
If Er. re-appears, press STOP and repeat the procedure.
If Er. continues to flash, replace control unit according to the instructions on pages 7 and 23. (State Er. number when requesting repairs).

2 → Proceed as described under 1, but without replacing the control unit, and consult after-sales service.






The heating muffle must be changed if heating performance becomes inadequate as a result of wear or if it is defective (see error indications ER23, E-24 on page 22).


Note:
In all procedures, avoid damaging the sealing ring (35), blade contacts (33) and switch pin (32).

1. Allow furnace to cool completely.
2. Press <key (furnace coast)
3. Once the furnace is fully closed, switch off 0/1 and disconnect mains plug.
4. Lift furnace heat (N) seat of the furnace base and place on a soft surface.
5. Remove heating muffle, following the instructions on page 25.
6. Mount new heating muffle, following the instructions on page 26.
7. Replace furnace heater furnace base, following the instructions on page 6.
8. Reconnect mains lead.


Removal of heating muffle




1. Press spring clip (36) inwards with a screwdriver and, using a second screwdriver, lift the end (36a) free. Remove spring clip.




2. Lever up and remove ring (37).




3. Lever out ends of insulation cord (38a) with spatula and gently pull out cord (38) by hand.




4. Carefully remove the stone lining segments (39).



5. Release the terminal screws (42a and 42b) of the terminals (41a and 41b) with approx. six turns of a screwdriver.



6. Withdraw end of wire (40b) from terminal (41b), lifting the heating muffle (40) slightly.




7. Withdraw end of wire (40a) from terminal (41a) and remove muffle (40).

Assembly of heating muffle

1. Clean inside wall of dome (30a) and recess in stone roof.

Note:


Handle wire ends and heating element with care. Avoid contact between heating element and fingers or hard objects.




2. Insert short end of wire (40a) into the terminal hole (41a).

Note:

Set mark on muffle opposite mark on interior wall of dome (30a).




3. Insert end of wire (40b) into terminal hole (41c). Set new heating muffle (40) in the recess in the stone roof.



4. Attach two ends of wire (40b) with terminal screws (42a) in the terminals (41a).
Note:
Both ends of the wire must be situated correctly and firmly secured in the terminals.

5. Reinsert the stone lining segments (39) so that they rest against the heating muffle (40). Introduce insulation cord (38) using a screw driver without tensioning along the entire circumference between the stone lining segments (39) and the interior wall of the cover (30a). Ensure that stone lining segments (39) fit tightly against the heating muffle. Lay ring (37) on the stone lining segments (39).




6. Secure spring clip (36) under the recess (30b) at the positions marked 'X'. Pull spring clip (36a) against the interior wall of the cover (30a) using two screwdrivers, as illustrated, until the ends (36a) and (36b) fit together. Release spring clip (36) with the two screwdrivers.
- Note:
The ends (36a and 36b) should now abut each other.

Furnace servicing

- Service is normally restricted to:
- occasionally checking that the steathed thermocouple is perpendicular
 - examining and cleaning the sealing ring (65) of the furnace base and the sealing ring (35) of the furnace head

Vacuum pump servicing

Refer to the operating instructions for the vacuum pump.



Summary of programs

Possible desired values in "Celsius mode"

*permanently programmed and stored (no display)
 **permanently programmed and stored (no display)

Program	Stand-by temperature	Temperature increase	Thawing temperature	Closing time	Holding time	Automatically set holding temperature	Automatically set holding time	Completion rate	Vacuumation	Vacuumation
P	°C	°C/min.	°C	minutes	minutes	°C	minutes	%	Yes/No	Yes/No
P01 to P08	Standard programs with defrost (automatically set) plus (cooling time 1 minute)									
P 01	125°	140°	40/1200	0.2-0.9 1-9	0.1-0.7 10-15			as P	1. Gamma closed**	and of H*
P 02	105°	140°	40/1200	0.3-0.9 9	hold = 0 without holding time (display 0)			up to 0 without completion control (no display)	without vacuum**	
P 03	105°	60°	40/1200	0.1-0.9 1-9	as P 01			as P 01	380°	start of H**
P 04	105°	300°	40/1200	0.1-0.9 9	as P 01			as P 01	without vacuum**	
P 05	125°	300°	40/1200	0.1-0.9 9	as P 01			as P 01	350°	start of H**
P 06	125°	300°	40/1200	0.1-0.9 1-9	as P 01			as P 01	without vacuum**	
P 07	105°	60°	40/1200	0.1-0.9 1-9	as P 01			as P 01	0-1 up to 0 without vacuum (no display)	Yes/No, 1 up to 1 start of H and of H display = 14 Vacuumation start of H
P 08	105°	300°	40/1200	0.1-0.9 1-9	as P 01			as P 01		
P 09 to P 15	Highly programmable programs (when Gamma 1 in mode)									
	100-700	30-140	150-200	0.1-0.9 1-9	as P 01			as P 01	as P 07/P 08	as P 07/P 08

Special programs with a second, automatically set holding temperature and holding time

P 20	100-700	30-40	300-1200	0.1-0.9	as P 01	125°	2	as P 01	as P 07/P 08	as P 07/P 08
P 21	100-700	30-40	300-1200	0.2-0.9	as P 01	150°	3	as P 01	as P 07/P 08	as P 07/P 08
P 22	100-700	30-40	300-1200	0.3-0.9	as P 01	180°	2	as P 01	as P 07/P 08	as P 07/P 08
P 23	100-700	30-40	300-1200	0.1-0.9	as P 01	100°	3	as P 01	as P 07/P 08	as P 07/P 08
P 24	100-700	30-40	300-1200	0.3-0.9	as P 01	125°	2	as P 01	as P 07/P 08	as P 07/P 08
P 25	100-700	30-40	300-1200	0.3-0.9	as P 01	125°	3	as P 01	as P 07/P 08	as P 07/P 08

Special program in which the (operating) cycles during time t_0 lasted with the first half with vacuum and in the second half without vacuum.

P76	100-200	30-140	300-1200	0.5-0.8	as P71			as P71	as P07 / P08	as P07 / P08
P77 or P87	Fixed programmatic special program (4 per furnace: 20 seconds)									
	200-300	30-140	300-1200	0.5-0.8 14.0	as P71			as P71	as P07 / P08	as P07 / P08
P88 or P90	Fixed programmatic night-time program (open furnace: 2 minutes; With 200-2020; 2 minutes; closes after opening at T = 450°C and does not continue to heat)									
	100-200	30-140	300-1200	0.3-0.9 2	as P71			as P71	as P07 / P08	as P07 / P08

Program for individual setting or completely switching off the boards:

P91	$V_1 = 0-1$ (use + sign);	$V_2 = 0-8$ (offset of signs); $t = 0-30$ sec for setting								
	Program for changing over from Celsius mode (°C) to Fahrenheit mode (°F) and vice versa.									
P92	For detailed description see page 13.									
P93	Program for silver test; for instructions see page 12.									
	403*	407*	405*							

Technical data

Power supply:

Single-phase AC
Standard model
Special models

Tolerated voltage fluctuations:

Power consumption:

Furnace and vacuum pump
Furnace alone (without pump)

Vacuum pump data (other manufacturers):

Max. permitted
rated current:
Max. permitted
peak current
(0.1 sec.)
Suction capacity
Final vacuum

220 V / 50-60 Hz
220 V / 50-60 Hz; 240 V / 50 Hz
118 V / 60 Hz
110 V / 60 Hz
+ 10% to -15%

approx. 1100 W (max. 1600 W)
max. 960 W

3 A at 200 to 240 V
4.4 A at 110 to 120 V
5 A at 200 to 240 V
8 A at 110 to 120 V

1.3 to 1.5 m³/h
27 to 40 mbar
(20 to 30 Torr)

Electrical fuses:

Values

200-240 V: T 6.3 A (heating circuit)
T 3.15 mA (mains)
T 3.15 A (pump)
110-118 V: T 12.5 A (heating circuit)
T 500 mA (mains)
T 5 A (pump)

Fuse dimensions

200-240 V: dia. 5 x 20 mm
110-118 V: dia. 6.3 x 32 mm

Dimensions of closed furnace:

width / depth / height

Effective firing chamber dimensions:

Max. firing temperature:

Weights:

Furnace head (K)

Control unit (S)

Furnace base (U)

Cooling plate (A)

Furnace complete

Space muffle

415 x 433 x 295 mm
dia. 80 mm, height 38 mm
1200°C (2192°F)

3.45 kg
2.30 kg
10.45 kg
0.30 kg
16.50 kg
0.25 kg

EC Declaration of Conformity

Document No./ Month, Year	105/12/95	Text:	a) Low Voltage Directive b) Electromagnetic Compatibility	Issued by:	Woclar AG FL-9494 Schaan
Manufacturer:	WOCCLAR Dental GesmbH Werk Burs Bremschstr. 16 A-6700 Bludenz-Burs		The appendices Form FN and Form BM-CP contain further information on the fulfillment of these Directives.	Place, Date:	Schaan, December 20, 1995
Address:				Valid Signature:	
Product Name:	Programat P95 Type Pt 1				
The product mentioned complies with the following European Directives:					
Number:	a) 73/23/EEC b) 89/336/EEC	The appendices are an integral part of this declaration. The declarator confirms the compliance with the Directives mentioned, but confers no warranty of attributes. The safety notes of the product documentation have to be observed.			

The Ivoclar Info System

This Info System comprises clearly arranged data and information for easy reference at your workplace.

Organisation of the System

- Foldable system binder with dividers.


Should you not yet possess an Ivoclar Info System, please contact any Ivoclar demonstrator or Ivoclar directly.



With the assortment you will receive a system binder containing a complete index of all parts. Each assortment is delivered with the respective punched Instructions for Use. Just remove the wire from the saddleback to file the Instructions for Use in the system binder.



Shorter Instructions for Use are perforated. After removal of the perforation, the Instructions for use can be filed in the system binder.



File Instructions for Use in the system binder.

This apparatus has been developed solely for use in dentistry. Scan up and operation should be carried out strictly according to the Operating Instructions. Liability cannot be accepted for damage resulting from misuse or failure to observe the instructions. The user is solely responsible for testing the apparatus for its suitability for any particular use as clearly stated in the instructions. Descriptions and data constitute no warranty of attributes and accuracy (including).

© 1997 Ivoclar Company Limited, Schaan
2212290/256/01/Av8/VD

Ivoclar Aktiengesellschaft
Bendererstrasse 2
FL-9494 Schaan / Liechtenstein
Tel. 075 / 235 35 35
Fax 075 / 235 33 60

Ivoclar- worldwide

Deutsche Ivoclar Dental GmbH
Postfach 1152
D-73471 Illwangen
Germany
Tel. 07961 / 889-0

Ivoclar Division France
B.P. 18
F-74410 Saint Jorioz
France
Tel. 50.88.64.00

Ivoclar-Vivadent s.r.l.
Via dell'Industria 16
I-39025 Naturno (BZ)
Italy
Tel. 0473 / 67.01.11

Ivoclar Comercial S.A.
c/Valdeiribas 82
E-28007 Madrid
Spain
Tel. 91 / 552 41 75

Ivoclar-Vivadent UK Limited
Meridian Soul
Leicester
LE2 2WY
Great Britain
Tel. 116 / 265 40 54

Ivoclar Dental (N.Z.) Ltd.
Wellesley Street
P.O. Box 5243
Auckland
New Zealand
Tel. (09) 670-5206

Ivoclar Pty. Ltd.
1 - 5 Overseas Drive
P.O. Box 307
Noble Park, Vic. 3174
Australia
Tel. 03 / 9295 95 90

Ivoclar Williams
Ivoclar North America, Inc.
23 Hannover Drive
St. Catharines, Ont. L2W 1A3
Canada
Tel. (800) 263-8182

Ivoclar Williams
Ivoclar North America, Inc.
175 Pineview Drive
Amherst, N.Y. 14228
USA
Tel. (800) 543-6825

Ivoclar-Vivadent
Branch Office in Poland
Swietokrzyska 36 apt. 37
PL-00-116 Warszawa
Poland
Tel. 22 / 620 24 64, 620 26 61

Dental Medico AB
Isafjordsgatan 9
Box 1194
S-16422 Kista
Sweden
Tel. 08 / 622 27 20

Ivoclar-Vivadent S.A. de C.V.
Av. Mazatlan No. 61, Piso 2
Col. Condesa
06170 México, D.F.
Mexico
Tel. (5) 553-0038 / 553-0972